

## Wenoll O<sub>2</sub> Circle Pack in carrying- and application bag



**image above:**  
Wenoll O<sub>2</sub> Circle Pack, closed, with carrying / shoulder strap

**left image:**  
Wenoll O<sub>2</sub> Circle Pack in use with a 1-litre Carbon cylinder, 300 bar (available on request)

Our approved and efficient Wenoll re-circulation system in a special carrying and application bag. Designed for connection to external oxygen sources with appropriate pressure regulators (such as CS 215) or to our new, extremely light 300-bar Wenoll systems WS 110, WS 120 and WS 150.

The Wenoll re-circulation system allows providing 100% oxygen concentration during longer periods - e.g. 7.5 hours with a 2 litre cylinder, 200 bar. The system is pre-assembled and, in a few simple steps, ready for service.

## ■ How to operate the O<sub>2</sub> Wenoll Circle Pack

- First remove the red cap from the Co<sub>2</sub> absorber and then connect the rebreathing bag through the opening of the bag from outside
- Remove the second red cap from the flex tube and connect the mask/mouthpiece
- Disconnect the O<sub>2</sub> connecting tube from the CO<sub>2</sub> absorber and connect it to the pressure regulator
- Follow the instructions in the Wenoll user manual to perform O<sub>2</sub> flush and -dosage

## ■ Contents of O<sub>2</sub> Wenoll Circle Pack

- 1 Co<sub>2</sub> absorber with O<sub>2</sub> connecting tube, 1,6m
- 1 tubing kit with overflow valve 2 mbar and flex tube
- 1 rebreathing bag, 3 litres
- 1 mask with airbag
- 1 mask holder for mask with airbag
- 1 single-use syringe for mask with airbag
- 1 silicone mouthpiece
- 1 nose clip
- 1 APL valve 0-60 mbar for artificial respiration
- 1 Wenoll system user manual
- 1 carrying and application bag with eyelets for fixing an O<sub>2</sub> cylinder



## ■ Dimensions and weights

External dimensions of carrying and application bag: \_\_\_\_\_ 32 x 27 x 13 cm (W x H x D)

Total weight: \_\_\_\_\_ 1.6 kg

## ■ Ordering information

Wenoll O<sub>2</sub> Circle Pack \_\_\_\_\_ Art.-Nr. 123 300

Thermo bag for rebreathing bag \_\_\_\_\_ Art.-Nr. 123 310